

## HPLC LIQUID FLOW METER - HMI

### Specification

- Desktop Architecture or in LAN network for your Labs.
- Provision for use within a SINGLE Location in a company.
- Each Lab can have one or two flow meters for HPLC Calibration process.
- User has a Login ID with Password protection.
- The HMI captures the readings on-line, displayed on LFM and transmitted through its RS 232 Serial COM port.
- HMI provides Date and time wise Report for all points of Calibration readings captured.
- The Report displays various flows, Actual flow rates, Mean value. There is date time stamp in the Header. Footer can have Analyst and Reviewer remarks.
- There is SET mode for programming company details like Company Name, Location, serial number, unit of measurement. The RUN mode is for normal operation.
- Operating License is provided for particular Serial number of the flow meter.
- The Unit can process Data as it is HMI with Mini PLC.
- The A4 size Formatted Report of readings can be printed on LAN printer if required.



MODEL : LFM 02H