

HIGH VOLUME CONCENTRATOR

High Volume Concentrator / N2 Evaporator Model : EV-PLUS-08

To accommodate 8 Nos. Beaker 200 ml capacity. Microprocessor based, water bath with adjustable temperature and timer for purging Nitrogen Gas with individually controlled sample concentration.

SS Rack accommodate Beakers as mentioned.

Rack Type	Working Volumes
50 x 70 mm - 100 ml Beaker	80 ml
60 x 80 mm - 150 ml Beaker	120 ml
70 x 95 mm - 250 ml Beaker	200 ml

PreCiVapTM
(PreCision Nitrogen Evaporator / Concentrator)



PreCi Analytic