



LCP 100

Specifications

- Flow Rate : 0.1 to 10.000 mL/min for 10 mL/min head
- Pressure : 0 to 6,000 psi
- Flow accuracy : $\pm 5\%$ for a flow rate of 0.1 to 10.000 mL/min, with 100% water @ 1000 psi
- Dimensions : 5.5" high x 10.375" wide x 17.5" deep
- Weight : 24 lb
- Power : 100-120 VAC, 50-60 Hz; or 220-240 VAC, 50-60 Hz
- Features : 4 line LCD display, Splash proof keypad, Prime purge valve, Auto stop timer, Service mode
- Optional : Pulse damper
- Nos of column : Upto 4 column can be washed simultaneously
- Remote inputs : RS-232