

Solid Sampling Hydraulic Press



Automatic Press
MODEL : HP-15TA



Manual Press
MODEL : HP-15TM



Mini-Pellet Press
MODEL : HP-Mini

KBr Die Set



Most commonly used Die for IR/FTIR for solid sampling of 13 mm pallet size consist on Anvil & Plunger, Top & Bottom Die Port, Extractor Ring, Oring.
Other sizes like 10 mm, 20 mm also available.

Dry Box (Thermostatic / PID Controller)



MODEL : DB-02



MODEL : DB-01

KBr Powder & Nujol Mull

The ideal method of storage of FTIR/XRF accessories which are affected by moisture since highly Hygroscopic in nature, which are to be stored in dry box.

- Input Voltage : 230 V AC, 50 Hz
- Dimensions : 407 L x 229 W x 368 H for **DB-02**
: 400 L x 210 W x 18 H for **DB-01**
- Temperature : Ambient +5°C to 70 °C-DB-02
Ambient +5°C to 45 °C-DB-01

Manual Press

- A Complete Laboratory hydraulic press producing a force about 15 tones use to make high quality 13 mm pallet used for IR / FTIR / XRF solid sampling.
- 15 ton laboratory hydraulic pallet press is a compact, elegant and robust machine, typically used by R&D & QC labs for various pelletizing applications for IR / XRF etc.
- The high pressure pumping unit supplies hydraulic fluid to the up-stroking ram of the cylinder. This causes the ram to rise steadily and positively in the upward direction. As a result, pressure is applied on any object placed between screw and piston top plate.

Mini-Pellet Press

- Full hydraulic operation and integrated pressure gauge
- Simple operation requiring minimal training
- Small and light weight
- Low cost of ownership
- Ideal for FTIR / KBr pellet preparation

Automatic Press

- Compact design and easy to operate electrical controls with pressing and ejection cycle.
- Pressing cycle has both auto and manual modes whereas ejection cycle has only manual (inching) mode
- Settable Parameters : Pressure, hold time & ram retracting time.
- Safety mechanisms to avoid over pressure & excessing ram stroke.
- Transperent polycarbonate safety guard has been provided in the area of operation
- Emergency Stop button to cut-off electrical supply (if required)
- Reputed MCB, overload relay & relay contactors are included in the Control-panel for electrical safety.
- Press is enclosed in a powder coated metallic cover, which is easily removable for maintenance purpose.



Pellet Holder



Agate Mortar and Pestle

The pellet holder is use to hold pallet (13 mm) of KBr, suitable to any IR / FTIR
Agate, Motar pestle use to prepare sample.

Liquid Sampling NaCl / KBr Windows



NaCl / KBr windows is available in circular & Rectangular shape.

Standard sizes :

Circular : dia 25 x 4 mm thickness

Rectangular : 38 x 19 x 4 mm

Other windows like AgCL & CaF₂ & other dimensions also available.

Universal Liquid Cell Holder (Demountable / Fixed Thickness)

Universal Liquid Cell Holder (Demountable / Fixed Thickness)



Universal Liquid Cell Holder



Mull Cell Holder

The cell holder is used to mount circular as well as rectangular windows & it can be used as Demountable Cell & Fixed Thickness Cell using different spacers supplied along with teflon washers of different sizes. Mull cell holder is used to mount circular window & for mull samples. Cell holder are supplied along with assorted spaces of size 0.1 mm, 0.2 mm, 0.5 mm & 1 mm.