

WELL PLATE CONCENTRATOR



MODEL : 96 X 2 WP

Nitrogen Concentrator/Evaporator 96 x 2 WELL PLATE

- Model : EV-PLUS 96 X 2 WP
- Capacity : One Standard 96 x 2 well plates Max.
- Sample volume : 1 mL (as per WP specs)
- Well plate Bay : 2 included as per requirement
- Number of Samples : Single Well plate of 96 samples
- Heating : Block Heating
- Timer Range : 1 to 99 minutes/indefinite Gas Regulator and Gas Gauge
- Range : 0 to 20 psi
- Gas Supply Requirements : Min. inlet pressure 60 psi, Max. inlet pressure 100 psi
- Exhaust Port : In-built ,
- Duct hose provided Power Requirements : 230 VAC at 50 Hz Net
- Weight : 15 kg Approx
- Technology : Vortex technology
- Gas control : Single input