



Model : EM UV-200



Printer

Technical Specifications

- Microcontroller based Electronic Module : For any instruments with Serial / USB / COM ports to print Data to Serial printers.
- Header generation -
 - (a. 35 characters per line including blanks)
 - (b. 2 lines programmable information)
 - (c. 3 lines fixed information)
- Display : 2 Lines LCD display, 16 Characters per Line.
- Footer generation : Automatic, 4 lines fixed information.
- Keypad : Feather touch Keypad.
- Modes of operation : Program, Set mode and Run mode.
- Ports for connection : 1 Input and 1 output port, D connector
- Calendar & Real time clock : Internal Day, Date DD-MM-YYYY and Time HH:MM:SS
- Communication protocol : Standard RS 232/485 serial Com port
- Baud Rate : 9600 / 2400 switch selectable.
- Power supply : 9V DC / 500 mA regulated supply.
- Programming Interface : Standard Windows HyperTerminal or compatible program.
- Data Count : Can set the readings from One to Max 10 readings .
- Data process Memory : Rounding off of 4th digit to match display.