



Gas Distribution Panel for GC, AAS & ICP

Gas Distribution Panel

- Gas distribution panel is used when there is no more than one gas chromatograph.
- It helps to control the flow & pressure of different gases for different GC's.
- Available in three different types. One is with pressure gauge, pressure regulator, ON/OFF valve. Second with pressure gauge and ON/OFF valve. Third is with only ON/OFF valve.



Gas Selector Box

Air N₂ H₂

Gas Selector Box

- Used for selection between two gases of different types.
- Available in two different types one with selector valve, pressure gauge & toggle valve, another with only two toggle valve.



Gas Distribution Panel Bench Model



Gas Distribution Panel Open Type

Gas Changeover Box (Cylinder vs Generator)

- Used for selection between cylinder & gas generator.
- To use gas from cylinder when gas generator is break down for service or maintenance.

PCI Analytical