




Installation Diagram for ICP




Gas Purification / Control System

- Gas purification panel is used to remove the vapours impurities like Moisture, CO, CO₂, Hydrocarbon from N₂ & Argon.
- Gas purification for AAS & ICP-MS is optional is one can also use gas distribution panel instead of gas purification panel
- CE certified
- With safety acrylic covers
- Inlet & outlet connection: 1/4" or 1/2"



Exhaust System for LCMS



- Blower 0.5 HP (pressurised), 3800 RPM
- SS Body, SS Duct
- Y shape suction connection
- SS Ball valve
- T Clamp
- MS bracket for blower