

◆ HPLC Solvent Waste Kit

Collect effluent from multiple sources while minimizing harmful solvent exhaust

- ◆ 10 port manifold - 1x1/2" NPT(F),3x1/4" NPT(F),and 6x1/4-28 UNF(F)-1 Each
- ◆ 83B Connection
- ◆ Polpropylene(PP) manifold
- ◆ PP Waste Collection Container , Carbouy 83mm Neck Size - 10 Ltr/20 Ltr
- ◆ SS Waste Collection Container,Carbouy 83mm Neck Size - 10 Ltr



pconlab[®]...
(Precision Consumables for Laboratory)