

Path Length	Capacity
10 mm	3.5 ml
1 mm	0.35 ml
2 mm	0.70 ml
5 mm	1.7 ml
20 mm	7.0 ml
30 mm	10.5 ml
40 mm	14 ml
50 mm	17.5 ml
100 mm	35.0 ml
10 mm (M) (Micro)	1.0 ml
10 mm (S) (Semi Micro)	1.4 ml
10 mm (U) (Ultra Micro)	0.7 ml
10 mm (U) (Ultra Micro)	0.5 ml
10 mm with Stopper Teflon with Round Bottom	3.5 ml

